

Miniaturization and Optimization of a Two-Plasmid Expression Pathway in myTXTL[®] Cell-Free Protein Expression Kit (Arbor Biosciences) Using Labcyte Echo[®] Acoustic Technology

Evelyn Eggenstein¹, Jared Bailey², Jean-Marie Rouillard¹, John Lesnick²

¹Arbor Biosciences, Ann Arbor MI, USA, ²Labcyte Inc. San Jose CA, USA

INTRODUCTION

In vitro transcription (TX)-translation (TL) is a rapidly developing technology with vast potential for synthetic biology and bio-engineering. myTXTL[®] (FIG. 1), an *E. coli* based cell-free protein expression system entirely relying on the endogenous TXTL machinery, is a versatile and all-in-one solution for many applications, such as the development of multi-stage gene circuits, and screening of complex protein libraries. For both examples, high sample throughput at an affordable cost is desirable.



FIG 1. myTXTL[®] workflow.

EXPERIMENTAL DESIGN

A typical myTXTL[®] reaction has a total volume of 12 μ L and consists of a ready-to-use myTXTL[®] Master Mix and a nucleotide template. Here, we pursued cell-free expression of a model protein (eGFP mutant) under transcriptional control of the T7 system. This requires the presence of T7 RNA polymerase, which was supplied to the myTXTL[®] reaction encoded on a plasmid (FIG 4).




FIG 4. *In vitro* protein production in myTXTL[®] using the T7 expression system.

RESULTS

Assay 1. Miniaturization of myTXTL[®] reaction

In the course of generating a standard curve for quantification of deGFP produced *in vitro*, the Echo[®] 525 Liquid Handler was able to accurately dispense a deGFP standard at different concentrations for volumes between 1 μ L and 12 μ L (FIG 5).




FIG 5. deGFP standard curve set up by acoustic transfer. deGFP fluorescence in the 384-well plate was acquired in a PHERAstar[®] reader at 485 nm (ex) and 520 nm (em).

Assay 2. Optimization of a two-plasmid myTXTL[®] reaction

Using acoustic liquid transfer, reagents were arrayed into 384-well plates to allow multiplexing of the dual-plasmid circuit (TABLE 1). Then, plates were incubated at 29 °C using a Labcyte MicroClima[®] lid to prevent evaporation, and fluorescence signal of deGFP produced in myTXTL[®] was continuously recorded over 16 hours.

TABLE 1. Experimental scheme for evaluating a two-plasmid myTXTL[®] reaction using Echo[®] 525 Liquid Handler for reagent transfer.

384-well plate matrix for multivariate DOE	P70a-T7rnap [pmol]											GFP STD Curve												
	100	75	50	37.5	25	10	7.5	5	2.5	1.25	0.62													
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
T7p14-deGFP [nM]	6.75	A-B																						
	5.63	C-D																						
	4.5	E-F																						
	3.38	G-H																						
	2.25	I-J																						
	1.13	K-L																						
	0.56	M-N																						
0.28	O-P																							
	Reagent											Acoustic Transfer Range												
	T7p14-deGFP											25-600 nL												
	P70a-T7rnap											25-400 nL												
	deGFP Standard											50-1000 nL												
	myTXTL [®] Master Mix											3000 nL												
	diH ₂ O											0-950 nL												

Maximum production of deGFP was observed in the presence of 4.5 nM T7p14-deGFP and 37.5 pmol P70a-T7rnap (FIG 6). At higher concentration of either plasmid, total deGFP yield was reduced. Alterations in the concentration of P70a-T7rnap were most pronounced over the time course of the experiment (FIG 7).




FIG 6. deGFP expression in myTXTL[®] as a function of plasmid concentration. Either T7p14-deGFP (blue) or P70a-T7rnap (purple) were kept constant, while the concentration of the corresponding plasmid was titrated over several magnitudes. Data recorded after 4 h (circle), 4.5 h (square), 5 h (triangle), 5.5 h (diamond) and 6 h (inverted triangle) is plotted.




FIG 7. Kinetic of deGFP production in myTXTL[®] under three different reaction conditions of the dual-plasmid system. The 4-6 h window is marked with dotted lines. Fluorescence signal was acquired every 5 min.

DISCUSSION

Arbor Biosciences' myTXTL[®] Master Mix combined with the Labcyte Echo[®] gives the ability to obtain fluorescent protein readout from a multivariable DNA input within a matter of hours. Rapid arraying, in combination with the reduced DNA and reagent demand, allows for a large increase in the number of possible experimental conditions as seen in Table 1. More work with acoustically arrayed myTXTL[®] will be performed to further highlight the flexibility of both systems in building a biological circuit.

HIGHLIGHTS

- Echo[®] 525 Liquid Handler reproducibly and reliably transfers myTXTL[®] Master Mix.
- Miniaturization reduces overall assay costs by a factor of three due to lowered reagent consumption.
- The speed of arraying allows for time-dependent reads of multiple plasmid inputs in a myTXTL[®] reaction.



FIG 2. Labcyte Echo[®] 525.

The Labcyte Echo[®] 525 Liquid Handler (FIG. 2) allows rapid, accurate and contact-free transfer of volumes at a nanoliter scale with acoustic sound. Its integrated Dynamic Fluid Analysis technology enables the Echo[®] system to easily adapt to various types of fluid without the need for calibration, which makes it ideal for handling multi-component solutions (FIG 3). This simplifies experimental setup and enables a maximum degree of flexibility for study design.




FIG 3. Acoustic droplet ejection (ADE) technology of the Echo[®] Liquid Handler.

Combining these two sophisticated platforms provides a comprehensive solution for fully-automated high-throughput *in vitro* protein expression: from template preparation to product analysis. In this study, we demonstrate the ability to process Arbor Biosciences' myTXTL[®] system in very low volumes and to accelerate DNA titration and dispensing using the Labcyte Echo[®] Liquid Handler. This method will lower overall cost by reducing assay time and reagent consumption, while increasing experimental complexity and minimizing the risk of contamination and human errors.