



Wheat Exome v1 Kit

Hybridization Capture for Targeted NGS

User Manual

Version 1.51 August 2021

Access additional manuals at arborbiosci.com/mybaits-manual

FOR RESEARCH USE ONLY. Not intended for diagnostic use.



© Daicel Arbor Biosciences
5840 Interface Drive, Suite 101
Ann Arbor, MI 48103 (USA)
+1 (734) 998-0751
info@arbor.daicel.com

MYBAITS is a registered trademark of Biodiscovery, LLC dba Arbor Biosciences. DAICEL ARBOR BIOSCIENCES is a registered trademark of Daicel Corporation. ION TORRENT is a registered trademark of Life Technologies Corporation. DYNABEADS is a registered trademark of Thermo Fisher Scientific. MYONE is a trademark of Thermo Fisher Scientific. NEXTERA, ILLUMINA, and TRUSEQ are registered trademarks of Illumina, Inc. KAPA is a registered trademark of Roche Molecular Systems, Inc. IDT and XGEN are registered trademarks of Integrated DNA Technologies, Inc. PACBIO is a registered trademark of Pacific Biosciences of California, Inc. OXFORD NANOPORE TECHNOLOGIES is a registered trademark of Oxford Nanopore Technologies, Ltd.



myBaits®

Wheat Exome v1 Kit Protocol

CONTENTS

INTRODUCTION

Procedure overview	3
Product compatibility	4
Changes since previous manual	4
Kit components and stability	4

REQUIREMENTS AND RECOMMENDATIONS

Input library	5
Blocking oligos	5
Equipment and reagents required	6

PROCEDURE Part 1: Hybridization Setup

W1.1 Prepare materials	7
W1.2 Hybridization Mix setup	8
W1.3 Blockers Mix setup	8
W1.4 Reaction assembly	9

PROCEDURE Part 2: Bind and Wash

W2.1 Prepare materials	10
W2.2 Wash Buffer X preparation	11
W2.3 Bead preparation	11
W2.4 Binding beads and hybrids	11
W2.5 Bead washing	11

PROCEDURE Part 3: Library Resuspension, Amplification, and Sequencing

W3.1 Prepare materials	12
W3.2 Enriched library recovery	12
W3.3 Library amplification	12
W3.4 Sequencing	13

APPENDIX


WA1 Troubleshooting	14
WA2 myBaits Wheat Exome v1 Protocol Quick Guide	15

INTRODUCTION

myBaits® is an in-solution NGS library target enrichment system, compatible with Illumina®, Ion Torrent®, and essentially any amplifiable sequencing library. We use a versatile nucleic acid synthesis technology to produce biotinylated RNA “baits” that are complementary to sequence targets, here roughly 200 Mbp of exon and exon-proximate space as defined by the version 1.1 annotation set of the IWGSC Chinese Spring reference genome sequence (version 1). Baits and other reagents for NGS target enrichment are supplied with your myBaits kit. After enrichment with myBaits, libraries may then be sequenced on the aforementioned platforms, or further prepared for PacBio® or Oxford Nanopore Technologies® sequencing.

Procedure overview

1. Amplified sequencing library, adapter blockers, and other hybridization reagents are combined.
2. Libraries are denatured, allowing blockers to hybridize to adapters. Baits are then introduced and hybridized to targets for several hours.
3. Bait-target hybrids are bound to streptavidin-coated magnetic beads and pulled out of suspension with a magnet.
4. Most non-target DNA is washed away. The remaining library is then amplified and either taken directly to sequencing or further treated.



Arbor Biosciences products compatible with this manual

This manual is compatible with Arbor Biosciences catalog numbers starting with 3091. Visit arborbiosci.com/mybaits-manual to download alternative and previous manuals, or email techsupport@arbor.daicel.com for assistance.

Changes since myBaits Wheat Exome v1 Manual version 0.9 (August 2019)

- **Block X** is now supplied and users are no longer required to supply their own blockers
- **Buffer E** is now supplied for post-cleanup bead resuspension

Kit components and stability

	Reagent	Cap color	Volume (8 rxn)	Volume (48+ rxn)
Box 1 Store at 4°C	Hyb N	Red	175 µL	500 µL
	Hyb S	Teal	750 µL	750 µL
	Beads	Clear	550 µL	1600 µL
	Binding Buffer	-	12 mL	36 mL
	Wash Buffer	-	20 mL	60 mL
Box 2 Store at -20°C	Hyb D	Yellow	70 µL	190 µL
	Hyb R	Purple	25 µL	70 µL
	Block O	Blue	90 µL	270 µL
	Block X	Orange	18 µL	106 µL
	Buffer E	Light Green	550 µL	1600 µL
Box 3 Store at -80°C	Baits	White	50 µL per 8 rxns	50 µL per 8 rxns

At the recommended storage temperatures, myBaits kit components have a shelf life of one year.

It is strongly recommended that sub-aliquots of Baits are made in reaction sizes appropriate for your experiment plans to minimize freeze-thaw cycles.

REQUIREMENTS AND RECOMMENDATIONS

Input library

TYPE Use myBaits with PCR-amplified and amplifiable Illumina TruSeq®-style or Nextera® Flex NGS libraries built from hexaploid, tetraploid, or diploid wheat. Inserts should be between 300 and 400 bp on average. Dual-indexed libraries are strongly recommended to reduce the hazard of mis-indexing induced by PCR jumping events. **It is NOT recommended to use myBaits with PCR-free libraries; additionally, myBaits are incompatible with libraries made using original Nextera or Nextera XT library preparation kits, or any library type containing biotin.**

AMOUNT Each Wheat Exome v1 enrichment reaction should be performed on a pool of 8 libraries (1 µg each; 8 µg total) concentrated to 7 µL with vacuum centrifugation or other means, and suspended in a nuclease-free buffer or water.

Blocking oligos

When ordering your myBaits kit, please indicate the sequencing library configuration you intend to enrich. The standard adapter blocking reagent provided with the kit (**Block X**) is compatible with Illumina TruSeq-style or Nextera-style libraries with single 6-12 bp or dual 6-12 bp indexing. For different adapter configurations, we recommend ordering Custom IDT® xGen® Blocking Oligos customized for your NGS library adapter sequences. At a concentration of 1 µg/µL, any custom adapter-blocking oligos can be used in lieu of myBaits Block X.

Equipment required

Item	Notes
50 mL tubes	1 per 44 reactions
Microcentrifuge tubes (1.5, 1.7, or 1.8 mL)	2 per 8 reactions
0.2 mL PCR strips with attached lids	2 per 8 reactions (e.g., VWR Cat# 93001-118)
Pipettors and tips for 0.5 - 500 μ L	Multichannel for 20 and 500 μL recommended
[†] Thermal cycler with heated lid compatible with 0.2 mL strips	1 or 2
Magnetic particle collector for microcentrifuge tubes	1 (e.g., ThermoFisher Cat# 12321D)
Magnetic particle collector for 0.2 mL strips	1 (e.g., Permagen Cat# S500)
Vortex mixer and mini-centrifuge for tubes and strips	
Water bath or incubation oven at 65°C	
Heat block for microcentrifuge tubes at 60°C	

[†] Ensure that the thermal cycler and strips allow no more than 4 μ L of 30 μ L volume evaporation overnight at 65°C

Reagents required

Reagent	Notes
Nuclease-free ("NF") water	900 μ L per reaction
PCR primers to amplify sequencing libraries after capture, e.g.:	
Illumina P5: AATGATACGGCGACCACCGA	2.5 μ L @ 10 μ M per reaction
Illumina P7: CAAGCAGAAGACGGCATACGA	2.5 μ L @ 10 μ M per reaction
PCR reagents for post-capture amplification	1 per reaction (e.g. Roche Cat# 07958927001)
PCR purification system, e.g., silica columns or SPRI beads	1 cleanup per reaction

PROCEDURE

PART 1: Hybridization setup

Pools of 8 sequencing libraries (1 µg each, for a total of 8 µg) are mixed with various blocking nucleic acids, denatured, and then combined with other hybridization reagents (including baits). These hybridization reactions incubate overnight to allow baits to encounter and hybridize with target library molecules.

W1.1 Prepare materials

Reagents	Equipment
Hyb N and Hyb S from Box 1	Nuclease-free microcentrifuge tubes (×2)
Hyb D and Hyb R from Box 2	0.2 mL strips with attached lids (×1 per 8 reactions)
Block O and Block X from Box 2	Pipettors and tips; multichannel for 20 µL recommended
Baits from Box 3 KEEP ON ICE	Vortex mixer and mini-centrifuge for above tube types
Libraries in pools of 8, 1 µg each, in 7 µL per reaction	Heat block set to 60°C
	Thermal cycler(s); 2 blocks recommended for >23 reactions

Program the thermal cycler:



Set lid temperature 5 to 10°C above each step temperature to minimize evaporation

Step	Temperature	Time
1	95°C	5m
2	65°C	5m
3	65°C	∞

W1.2 Hybridization Mix setup

- Once the Hyb reagents have thawed, vortex to homogenize and then briefly centrifuge.



Heat Hyb N and Hyb S to 60°C and vortex to dissolve any precipitate present after thawing

- Assemble the Hybridization Mix in a microcentrifuge (MC) tube, briefly vortex and briefly centrifuge the contents to collect. **The following volumes are already adjusted for pipetting error:**

Component	μL / Reaction
Hyb N	9.25
Hyb D	3.5
Hyb S*	0.5
Hyb R	1.25
Baits	5.5
TOTAL	20

**Cloudiness caused by Hyb S addition will clear after step 3*

- Incubate the Hybridization Mix at 60°C for 10 minutes in the heat block. Vortex occasionally to collect the condensate. Remove from the heat block and let sit 5 minutes at room temperature.
- For each capture reaction, aliquot 18.5 μL of Hybridization Mix to a 0.2 mL well/tube.

These reaction aliquots of Hybridization Mix are now referred to as "HYBs"

W1.3 Blockers Mix setup

- Assemble the Blockers Mix in an appropriately-sized tube and mix by pipetting. **The following volumes are already adjusted for pipetting error:**

Component	μL / Reaction
Block O	5.4
Block X	2.1
NF Water	-
TOTAL	7.5

- For each capture reaction, aliquot 7 μL of Blockers Mix to a 0.2 mL well/tube.
- Add 7 μL of individual or pooled libraries to each Blockers Mix aliquot and mix by pipetting.

These libraries mixed with Blockers Mix aliquots are now referred to as "LIBs"

W1.4 Reaction assembly


Double-check the thermal program:




Set lid temperature 5 to 10°C above each step temperature to minimize evaporation

Step	Temperature	Time
1	95°C	5m
2	65°C	5m
3	65°C	∞


- Put the **LIBs** in the thermal cycler, close the lid, and start the thermal program.



- Once the cycler reaches 65°C during step 2, pause the program, put the **HYBs** in the thermal cycler, close the lid, and resume the program.



- After step 2 of the program is complete, leaving all tubes in the thermal cycler, **pipette 18 µL of each HYB to each LIB**. Use a multichannel pipettor for easier execution. Gently homogenize by pipetting up and down 5 times.



- Dispose of the HYB tubes. Briefly spin down the LIBs, return to the thermal cycler, close the lid, and allow the reactions to incubate overnight (16 to 24 hours).

PART 2: Bind and Wash (“Cleanup”)

After overnight incubation, bait-target hybrids are bound to streptavidin-coated magnetic beads, and then most non-target DNA is removed with several rounds of washing with a warm buffer.

W2.1 Prepare materials

 *Start at least 90 minutes before intended hybridization stop time*

Reagents

Hyb S (Box 1) *

Binding Buffer (Box 1) *

Wash Buffer (Box 1) *

Beads (Box 1)

Nuclease-free (NF) Water (up to 900 µL per reaction)

** Allow these reagents to come to room temperature before use; warm to 60°C and vortex to dissolve precipitate if necessary*

Equipment

Water bath or incubation oven set to 65°C

Receptacles for 50 mL tubes, 0.2 mL strips and microcentrifuge tubes compatible with above incubation device

Vortex mixer and mini-centrifuge for 0.2 mL strips and MC tubes

Magnetic particle collector(s) (MPC) for above strips and/or tubes

When using only a microcentrifuge (MC) tube-compatible MPC	When using a 0.2 mL tube-compatible MPC
Nuclease-free 50 mL tube, 1 per 44 cleanups	Nuclease-free 50 mL tube, 1 per 68 cleanups
Nuclease-free MC tubes, 1 per reaction	Nuclease-free 0.2 mL PCR strips with individually-attached lids, 1 vessel per reaction
Heat block set to 65°C	Thermal cycler set to 65°C
Pipettors and tips for 20 – 500 µL	Pipettors and tips for 20 – 200 µL; multichannel pipettor strongly recommended

W2.2 Wash Buffer X preparation

This step generates enough Wash Buffer X for 44 reactions in microcentrifuge (“MC”) tube cleanup format, and 68 reactions in 0.2 mL cleanup format; scale up or down if needed.

1. Thaw and thoroughly homogenize Wash Buffer and Hyb S prior to aliquoting in order to dissolve any visible precipitate; warm slightly if necessary.
2. Combine 400 μ L Hyb S, 39.6 mL NF water and 10 mL Wash Buffer in a 50 mL tube. Vortex thoroughly, label “Wash Buffer X.” *Wash Buffer X can be stored at 4°C for 1 month*
3. Heat the Wash Buffer X to 65°C in the water bath or oven for at least 30 minutes before use.

W2.3 Bead preparation Prepare beads immediately prior to use

1. For each capture reaction, aliquot 30 μ L beads to a microcentrifuge tube.
2. Pellet the beads in the MPC until the suspension is clear (1-2 minutes).
Leaving the tubes on the magnet, remove and discard the supernatant.
3. Add 200 μ L Binding Buffer to each bead aliquot. Vortex to resuspend the beads and centrifuge briefly. Place tube in the MPC and pellet beads; remove and discard the supernatant.
4. Repeat Step 3 above twice for a total of three washes.
5. Resuspend each washed bead aliquot in 70 μ L Binding Buffer. If proceeding to washing in 0.2 mL format, transfer the aliquots to PCR strips.

TIP: Beads can be prepared in 8 (or fewer) reaction batches (240 μ L) in a microcentrifuge tube. Multiply all volumes by the number of reactions in the batch; i.e., for 8 reactions-worth, wash with 1.6 mL and resuspend in 560 μ L Binding Buffer, then aliquot 70 μ L suspension to individual tubes.

W2.4 Binding beads and hybrids

1. Heat the bead aliquots to 65°C for at least 2 minutes.
2. Transfer each capture reaction to the heated bead aliquots. Mix by pipetting.
3. Incubate the libraries + beads on the heat block or thermal cycler for 5 minutes. Agitate at the 2.5 minute mark by flicking or inverting the tubes to keep the beads suspended, followed by briefly centrifuging.

W2.5 Bead washing

1. Pellet the beads with the MPC until the solution is clear. Remove and discard the supernatant.
2. Add 375 μ L (MC tube format) or 180 μ L (0.2 mL format) warmed Wash Buffer X to the beads, remove from the MPC, place on heat block for 15 seconds, and briefly vortex or mix by pipetting. Briefly centrifuge the mixture.
3. Incubate for 5 minutes at 65°C in the heat block or thermal cycler. Agitate at the 2.5 minute mark via gentle vortexing and then briefly centrifuge.
4. Repeat steps 1 through 3 two times for MC tube format (three washes total), or three times for 0.2 mL format (four washes total). **After the last wash and pelleting, remove as much liquid as possible without touching the bead pellet.**

PART 3: Library Resuspension, Amplification, and Sequencing

Bead-bound enriched library is resuspended in Buffer E and amplified.

W3.1 Prepare materials

Reagents	Equipment
Buffer E (Box 2)	Tubes appropriate for PCR master mix assembly
PCR primers for amplifying libraries (e.g., P5 and P7)	Tubes or strips for 50 µL PCR amplification
PCR reagents for post-capture amplification	Pipettors and tips capable of 5 – 100 µL volumes
PCR purification system, e.g., silica columns or SPRI beads	Vortex mixer and mini-centrifuge for above tube types
Nuclease-free (NF) Water	Thermal cycler

W3.2 Enriched library recovery

1. Add 30 µL Buffer E to the washed beads and thoroughly resuspend by pipetting.
Then, depending on your library amplification system, choose workflow A or B:

WORKFLOW A: When using KAPA HiFi HotStart or NEB Ultra II Q5 polymerase systems for amplification

- 2A. Proceed directly to section S3.3 using this bead resuspension as the template in amplification

WORKFLOW B: When not using the polymerase systems for amplification in workflow A

- 2B. Incubate the suspension at 95°C for 5 minutes
- 3B. Immediately pellet the beads in the MPC and collect the supernatant containing the enriched libraries

W3.3 Library amplification

This is an example post-capture amplification using KAPA HiFi HotStart ReadyMix and Illumina libraries:

1. Assemble the following PCR master mix:

Component	Final Concentration	µL / Reaction	Sequence
NF Water	-	5	-
2X KAPA HiFi HotStart Ready Mix	1 X	25	-
P5 library primer (at 10 µM)	500 nM	2.5	AATGATACGGCGACCACCGA
P7 library primer (at 10 µM)	500 nM	2.5	CAAGCAGAAGACGGCATACGA
Enriched Library (on- or off-bead)	-	15 *	-
TOTAL		50	-

*Remaining bead-bound library can be stored at -20°C for several months.

W3.3 Library amplification (continued)

- Cycle the reactions with the following thermal program:

Step	Temperature	Time	
1	98°C	2 minutes	
2	98°C	20 seconds	×10 cycles
3	60°C	30 seconds	
4	72°C	45 seconds	
5	72°C	5 minutes	
6	8°C	∞	

- After amplification:
 - **If beads were included in the amplification reaction and you intend to use silica columns for purification**, pellet the beads first and purify only the supernatant.
 - Otherwise, purify the reaction using your preferred PCR cleanup (e.g., silica columns or SPRI beads).

The enriched libraries are now ready for quantification, quality-assessment, and sequencing. Post-capture yields are typically between 50 and 300 ng as assessed with quantitative PCR.

W3.4 Sequencing

To sequence the targets in hexaploid wheat to 30X unique average unique read depth, we recommend 18 Gbp (60 M read-pairs) PE150 sequencing per library (144 Gbp or 480 M read-pairs per pool of 8 libraries). Tetraploid and diploid wheat require 12 or 6 Gbp PE150 sequencing, respectively.

APPENDIX

WA1: Troubleshooting

My on-target read percentage (before PCR duplicate collapse) is lower than 60%

Ensure that 8 micrograms total library material (as measured with quantitative PCR or intercalating dye such as) has been used per enrichment reaction. Lower input amounts and/or use of other adapter-blocking oligos than Block X may reduce percent on-target significantly.

During hybridization, my thermal cycler dropped below 65°C

This kit should return between 50 and 70% reads on-target when enriching libraries made from hexaploid wheat. If the temperature drops during hybridization, you can expect a lower on-target read proportion and target read complexity, but not outright enrichment failure. Shallow preliminary sequencing will determine whether targets are likely to be retrieved at sufficient coverage within budget.

I observe a high ratio of PCR duplicates in my enriched library sequence data

Sequencing 18 Gbp per library using PE150 protocol (i.e., 60 M read-pairs per library) on Illumina platforms typically yields fewer than 30% PCR duplicates following alignment to the IWGSC Chinese Spring reference genome assembly. If you observe higher PCR duplicate rate at 18 Gbp (or lower) sequencing depth, ensure 1 ug of library as measured with verified library quantitative PCR was used per enrichment reaction. Also ensure the libraries used the manufacturer's recommended nanograms input DNA per prep.

WA2: myBaits Procedure Quick Guide - Wheat Exome V1

- For each reaction, build the following Mixes; pipetting error is built in:

Hybridization Mix	
Component	μL / Reaction
Hyb N	9.25
Hyb D	3.5
Hyb S	0.5
Hyb R	1.25
Baits	5.5
TOTAL	20

Blockers Mix	
Component	μL / Reaction
Block X	2.1
Block O	5.4
NF Water	0 ⁺
TOTAL	7.5

- After pre-warming the Hybridization mix for 10 minutes @ 60°C, for each reaction, aliquot **18.5 μL** of Hybridization Mix to their own tubes – now **“HYBs”**.
- For each reaction, aliquot **7 μL** of Blockers Mix and then add **7 μL** of each library – now **“LIBs”**.
- Incubate the LIBs in the thermal cycler for 5 minutes @ 95°C and then drop to 65°C. **Be sure to use a heated lid set to 5 to 10 °C higher than each step temperature.**
- Put the HYBs in the thermal cycler and warm to 65°C for 5 minutes.
- Transfer **18 μL** of each HYB to each LIB, mix by pipetting, and incubate overnight.
- 1.5 hours before step 9**, prepare Wash Buffer X by combining 400 μL Hyb S, 39.6 mL nuclease-free molecular biology-grade water and 10 mL Wash Buffer in a 50 mL tube. Vortex thoroughly and warm to 65°C for at least 45 minutes.
- Prepare 30 μL of beads per reaction by washing three times in 200 μL Binding Buffer. Resuspend the washed bead aliquots in 70 μL Binding Buffer and warm the suspensions to 65°C for at least 2 minutes.
- Combine the warmed beads with the hybridization reactions and incubate for 5 minutes at 65°C, agitating at 2.5 minutes to keep beads suspended.
- Pellet the beads and remove the supernatant. If using microcentrifuge tubes for cleanup, wash the beads three times with 375 μL warmed Wash Buffer X, incubating 5 minutes at 65°C. Wash four times with 180 μL washes if using a 96-well magnetic particle concentrator and 0.2 mL strips/tubes.
- Resuspend the beads in 30 μL Buffer E and then use 15 μL of this in a 50 μL amplification reaction with KAPA[®] HiFi or NEB Ultra II Q5 polymerase systems. If not using these polymerase systems, instead elute the library from the beads by incubating the suspension for 5 minutes at 95°C, immediately pellet the beads, and then use 15 μL of the supernatant in a 50 μL amplification reaction.
- Purify the amplification reactions using silica columns or SPRI beads. **If using silica columns and beads were included in the amplification reaction, pellet the beads first and purify only the supernatant.**