# Synthetic Oligonucleotide in situ Hybridization Probes **Provide Flexible, High-Fidelity Target Detection**

### Linda Barthel, M.Sc., Alison Devault, Ph.D., Tianying Lan, Ph.D., Benjamin Steil, Ph.D. | Daicel Arbor Biosciences

The use of in situ hybridization is a long-standing application used to visualize a wide range of DNA or RNA sequences in chromosomes, cells, tissue sections or whole tissue samples. There are several technical challenges and limitations associated with traditional probe sources, such as BACs and other biologically-derived material. Many of these problems are overcome by using oligonucleotides synthesized in a pool as a probe source. myTags in situ hybridization probes from Daicel Arbor Biosciences are bioinformatically-designed, fully customizable labeled or immortalized oligo pools ready for in situ hybridization in virtually any species or research model organism. Highlighted are several published and ongoing applications using myTags, including in visualization of genome organization and chromatin topology, chromosome identification, scaffold assembly, haplotyping, detection of large-scale structural variation, target gene mapping and expression detection. The design, synthesis, and downstream treatment processes have been streamlined to maximize specificity and flexibility, while minimizing cost.

## Key Applications

Visualizing target sequences on chromosomes, in cells and nuclei by fluorescence in situ hybridization (FISH) is a fundamental tool in genomics and cell biology. Recently, in situ hybridization applications have benefited tremendously from synthetic oligo pools carefully designed to maximize specificity and sensitivity in a variety of hybridization experiments. myTags in situ probes are a highly flexible and cost-effective toolkit deploying this transformative technology. Highlighted here are numerous examples of their successful use in the literature:

Applications	Description	Featured Publications
3D Visualization and Analysis of Genome Organization	Coupling oligo-FISH with 3D microscopy and image reconstruction to study spatial organization of the nuclear genome, chromatin topology and interactions between chromosomes/chromatin.	Giorgetti, L. <i>et al.</i> (2016), <i>Nature</i> , doi: <u>10.1038/nature18589</u> . Zheng, M. <i>et al.</i> (2019), <i>Nature</i> ,doi: <u>10.1038/s41586-019-0949-1</u> . Boyle, S. <i>et al.</i> (2020), <i>Genes &amp;</i> <i>Development</i> , doi: <u>10.1101/gad.336487.120</u> .
Super-resolution Imaging of Nuclear Organization and Structures	Oligo-FISH in combination with super-resolution imaging techniques allows for analyzing the subcellular localization and spatial arrangement of targeted DNA sequences and RNA transcripts at the single cell level.	Beliveau, B.J. et al. (2015), Nature Communications, doi: <u>10.1038/ncomms8147</u> . Wang, P. et al. (2020), Genome Biology, 21(1), doi: <u>10.1186/s13059-020-02030-2</u> Kubalová, I. et al. (2021) , Cell Biology, doi: <u>10.1101/2021.09.16.460607</u> .
Chromosome Identification or Indexing, Karyotyping and Chromosome Painting	Development of chromosome identification/indexing system with oligo-FISH allows for karyotyping based on individually identified chromosomes and studies of chromosome-scale genetic adaptation and evolution. Chromosome painting using chromosome or haplotype-specific oligo probes for studies of chromosomal evolution, monitoring chromosome pairing and transmission and assessing genomic sequence assembly.	do Vale Martins, L. et al. (2019), Nature Communications, doi: <u>10.1038/s41467-019-12646-z</u> . Agrawal, N. et al. (2020), Frontiers in Plan Science, doi: <u>10.3389/fpls.2020.598039</u> . de Oliveira Bustamante, F. et al. (2021), Theoretical and Applied Genetics, doi: <u>10.1007/s00122-021-03921-z</u> .
Detection of Chromosome Structure and Rearrangement	FISH-based assays to reveal chromosome structures and rearrangements including translocations, segmental duplications/deletions, and presence-absence variations, which are essential for studying recombination between chromosomes, chromosome synteny and evolution, massive gene duplication and amplification events, etc.	Albert, P.S. et al. (2019), PNAS, doi: <u>10.1073/pnas.1813957116</u> . Piperidis, N. and D'Hont, A. (2020),, The Plant Journal, doi: <u>10.1111/tpj.14881</u> Wang, K. et al. (2022), Chromosome Research, doi: <u>10.1007/s10577-021-09680-3</u> .
Evaluation of Gene Expression and Regulation	Visualization of chromatin interactions and spatial arrangement/localization of target genes, nascent or mRNA transcripts. Using oligo-FISH allows for investigation of transcriptional and epigenetic regulation activities including the detection of SARS-CoV-2 RNA.	Wang, P. et al. (2020), Genome Biology, doi: <u>10.1186/s13059-020-02030-2</u> . Markaki, Yolanda and Plath, Kathrin (2021), doi:10.1016/j.cell.2021.10.022. Mirabelli, C. et al. (2021), <i>PNAS</i> , doi: <u>10.1073/pnas.2105815118</u> .

Coming soon! Catalog whole chromosome paint probes for human & mouse





phone: 1-734-998-0751 twitter: @ArborBio



my Tags