

myTXTL Cell-Free Expression of Antibodies and Antibody Fragments

Nathan Gardner, Jacob Krol, Jacob Enk, Ben Steil

Daicel Arbor Biosciences, Ann Arbor, MI, USA



Abstract

Antibody discovery and screening processes are time, labor, and cost intensive primarily due to challenging protein expression. A novel myTXTL[®] Antibody/DS Master Mix enables rapid, cell-free expression of antibody constructs to greatly reduce overall time and resources required for antibody discovery. Historically, cell-free protein expression systems suffered from low yields, inconsistent batch quality and an inability to synthesize active antibody constructs due to a reducing environment. We demonstrate the lot-to-lot consistency of the myTXTL Antibody/DS Master Mix for the expression of a functional disulfide bond containing enzyme. Then, we demonstrate that myTXTL Antibody/DS Master Mix outperforms several existing cell-free systems on the market for both IgG and Fab expression. Due to the cell-free environment of myTXTL, dilutions of endpoint reactions can be used directly in ELISA assays for rapid antibody evaluation. An ELISA of Trastuzumab IgG expressed using several cell-free systems revealed that myTXTL Antibody/DS Master Mix produced 3-fold more active IgG than the nearest competitor. Omitting *in vivo* expression and consistently delivering industry-leading yields enables myTXTL Antibody/DS Master Mix to accelerate antibody screening and evaluation workflows.

Methods

GLucDura, Trastuzumab IgG and Fab fragment were expressed from plasmids containing a T7 promoter. Target proteins were purified using either MagStrep beads (IBA Lifesciences) or Protein A beads (GenScript). Trastuzumab IgG ELISA was conducted with an ELISA kit from IBL America (Cat. No. TM09013). Cell-free reaction parameters are summarized in Table 1. myTXTL reactions contained target protein plasmid(s), a helper plasmid expressing T7 RNA polymerase, and myTXTL Antibody/DS Master Mix. Competitor kit reactions were setup according to manufacturer's recommendations. Reaction setup and incubation conditions are summarized in Table 1.

Table 1: Cell-free reaction setup

Cell-Free System	# Components per reaction	Reaction Volume (uL)	Reaction Time (hours)	Reaction Temperature (°C)	Shaking Speed (RPM)
myTXTL	3	12	16	27	0
G	9	12	24	37	0
N	7	50	24	25	1000
L	2	50	48	25	700

Lot-to-Lot Consistency of myTXTL Antibody/DS Master Mix





Figure 1. GLucDura expression with Antibody/DS Master Mix: Three independent lots of myTXTL Antibody/DS Master Mix were evaluated for expression of active *Gaussia*-Dura luciferase (GLucDura) using a luminescent activity assay (NanoLight). Expression was tested in two independent rounds.

- All 3 lots of Antibody/DS Master Mix yielded significant and consistent amounts of active GLucDura enzyme, a protein with 5 disulfide bonds.

Conclusions

- myTXTL Cell-Free Expression Antibody/DS Master Mix has **robust lot-to-lot performance** for expression of proteins bearing disulfide bonds.
- The Antibody/DS Master Mix expressed **high yields of functional Fab and IgG**, outperforming competitors.
- Use of the myTXTL cell-free expression system can reduce the cost and time of antibody discovery workflows by enabling **antibody expression and evaluation within 24 hours**.
- Want to learn more about cell-free protein expression? Visit our **Booth #711!**

myTXTL Cell-Free Expression Antibody/DS Workflow



Cell-Free Expression and Evaluation of Trastuzumab Fab

Trastuzumab Fab was expressed in myTXTL Antibody/DS Master Mix and competitor kits with light and heavy chains under a single operon. Expressed Fab was incubated with HER2 antigen and pulled down by a Twin-Strep-tag[®] on the heavy chain to demonstrate antigen binding.




Figure 2. Evaluating cell-free expression of Trastuzumab Fab: SDS-PAGE visualization of pull-down eluates under non-reducing and reducing conditions reveals antigen and native Fab in all cell-free systems tested (Fig. 2A). Concentrations of Fab were determined from SDS-PAGE bands using a BSA loading control (Fig. 2B).

- myTXTL expressed high yields of Fab, more than competitors N and G (3-fold and 6-fold) and pulled down more HER2 antigen
- Bead purified Fab was of high purity for all kits tested

Twin-Strep-Tag[®] is a registered trademark of IBA Lifesciences

ELISA and Protein A Purification of Cell-Free Expressed Trastuzumab IgG

Trastuzumab IgG was expressed from two separate plasmids at a 1:1 ratio (HC:LC) in myTXTL Antibody/DS Master Mix for 24 hours and in 3 competitor kits. Duplicate myTXTL and competitor endpoint reactions were diluted and used immediately for an anti-HER2 ELISA.




Figure 3. Evaluating cell-free expression of Trastuzumab IgG. Four dilutions of each reaction were evaluated by ELISA to determine average IgG concentrations in endpoint reactions (Fig. 3A). Protein A magnetic bead purification of Trastuzumab IgG expressed using myTXTL Antibody/DS Master Mix shows that the 154 kDa IgG band is an abundant product with smaller species also evident on the gel (Fig 3B).

- myTXTL expressed high yields of functional IgG, 3x to 11x more IgG per unit volume than competitors
- IgG expressed in myTXTL can be purified with Protein A beads

myTXTL Antibody/DS Master Mix Features:

Template DNA
Plasmid or Linear

Reaction Volume
Scalable

Antibody Yields
Industry-leading

