

# myTXTL<sup>®</sup> Protocol

## His-tag Purification of Proteins from myTXTL<sup>®</sup> Reactions

**Beads:** HisPur™ Ni-NTA Magnetic Beads (Cat. No 88831, ThermoFisher)

**Binding Buffer** Composition:

- PBS (Phosphate buffered saline)

**Wash Buffer 1** Composition:

- 100 mM Na<sub>2</sub>HPO<sub>4</sub>-NaH<sub>2</sub>PO<sub>4</sub> (pH 8.0) (or potassium phosphate equivalent)
- 450 mM NaCl

**Wash Buffer 2** Composition:

- 100 mM Na<sub>2</sub>HPO<sub>4</sub>-NaH<sub>2</sub>PO<sub>4</sub> (pH 8.0) (or potassium phosphate equivalent)
- 600 mM NaCl
- 10-20 mM Imidazole (pH 8.0)

**Elution Buffer** Composition:

- 100 mM K<sub>2</sub>HPO<sub>4</sub>-KH<sub>2</sub>PO<sub>4</sub> (pH 8.0) (or sodium phosphate equivalent)
- 600 mM NaCl
- 350 mM Imidazole (pH 8.0)

### His-Bead Purification Protocol

|                           |   |
|---------------------------|---|
| <b>Sample Preparation</b> | <ol style="list-style-type: none"><li>1. Following myTXTL expression, centrifuge myTXTL reactions at 16,000 xg for 3 minutes to pellet any insoluble material. Then, place reactions on ice.<ol style="list-style-type: none"><li>a. If insoluble material is a significant concern, transfer the supernatant into a new tube</li></ol></li><li>2. Dilute 10-20 µL of myTXTL reaction product to 100 µL with Binding Buffer and place at room temperature (RT).</li></ol>   |
| <b>Bead Equilibration</b> | <ol style="list-style-type: none"><li>1. For Pro Master Mix samples, transfer 40 µL of the stock Bead slurry into a 1.5-mL tube and for Antibody/DS Master Mix samples, transfer 20 µL.<ol style="list-style-type: none"><li>a. Ensure the stock Bead slurry is mixed well before transfer.</li></ol></li><li>2. Add 160 µL Binding Buffer to the Bead slurry and vortex for 10 seconds at mid power.</li><li>3. Pulldown the beads with a magnetic stand and discard the supernatant.</li><li>4. Add 200 µL Binding Buffer to the Beads and vortex the beads for 10 seconds at mid power.</li><li>5. Pulldown the beads with a magnetic stand and discard the supernatant.</li></ol> |
| <b>Sample Binding</b>     | <ol style="list-style-type: none"><li>1. Add 100 µL Sample to the beads and vortex for 5 seconds at medium-low power.</li><li>2. Shake the sample with beads for 10 minutes in a ThermoMixer at 1400 RPM, RT.<ol style="list-style-type: none"><li>a. Alternately, tap the sample periodically to suspend beads.</li></ol></li><li>3. Pulldown the beads with a magnetic stand and discard the supernatant.</li></ol>   |
| <b>Washing</b>            | <ol style="list-style-type: none"><li>1. Add 400 µL Wash Buffer 1 to the Bead sample and vortex for 10 seconds at mid power.</li><li>2. Pulldown the beads with a magnetic stand and discard the supernatant.</li><li>3. Add 400 µL Wash Buffer 2 to the Bead sample and vortex for 10 seconds at mid power.</li><li>4. Pulldown the beads with a magnetic stand and discard the supernatant.</li><li>5. Add 200 µL Wash Buffer 2 to the Bead sample and vortex for 10 seconds at medium-low power.</li><li>6. Pulldown the beads with a magnetic stand and discard the supernatant.</li></ol>  |

## His-Bead Purification Protocol *cont.*

|                |  |
|----------------|--|
| <b>Elution</b> | <ol style="list-style-type: none"><li>1. Add 25 <math>\mu</math>L Elution Buffer to the Bead sample and vortex briefly at mid power.</li><li>2. Shake the sample with beads for 5 minutes in a ThermoMixer at 1450 RPM, RT.<ol style="list-style-type: none"><li>a. Alternately, tap the sample periodically to suspend beads.</li></ol></li><li>3. Pulldown the beads with a magnetic stand and transfer the eluate to a new 1.5-mL tube.<ol style="list-style-type: none"><li>a. This eluted volume should contain your His-tagged protein.</li><li>b. Load 2-9 <math>\mu</math>L of the eluate in an SDS-PAGE gel to visualize your target protein.</li></ol></li></ol> |
|----------------|--|

## Appendix: Endotoxin Removal

Endotoxin levels can be severely depleted, to about 10-20 EU/mL, using a simple wash procedure during HisPur Magbead Purification

|                                  |  |
|----------------------------------|--|
| <b>Materials</b>                 | <p><b>Detergent:</b> TritonX-114</p> <p><b>Triton Wash Buffer</b> Composition:</p> <ul style="list-style-type: none"><li>- 100 mM Na<sub>2</sub>HPO<sub>4</sub>-NaH<sub>2</sub>PO<sub>4</sub> (pH 8.0) (or potassium phosphate)</li><li>- 450 mM NaCl</li><li>- 1.25% Triton (%v/v)</li></ul> <p>Store at 4°C</p>  |
| <b>Endotoxin Removal Washing</b> | <ol style="list-style-type: none"><li>1. During HisPur purification, replace the first wash buffer with 400 <math>\mu</math>L Triton Wash Buffer and incubate the sample at 4-10°C for 10 minutes using a ThermoMixer at 1400 RPM.</li><li>2. Pulldown the beads with a magnetic stand and discard the supernatant.</li><li>3. Perform the standard HisPur washing steps for Wash 2 and 3.</li></ol> |

For research use only. Daicel and Daicel Arbor Biosciences are registered trademarks of Daicel Corporation.

Questions? Contact us via the methods listed below. Our team is happy to assist you!



**Web:** [www.arborbiosci.com](http://www.arborbiosci.com)  
**Email:** [info@arbor.daicel.com](mailto:info@arbor.daicel.com)  
**Phone:** 1-734-998-0751  
**X:** @ArborBio  
**Bluesky:** @arborbio.bsky.social



Rev. Date 02/2026