PalaeoChip A capture enrichment approach to ancient environmental DNA

Tyler J. Murchie¹, Jennifer Klunk², Melanie Kuch¹, Ana Duggan¹, Marissa L. Ledger⁴, Emil Karpinski¹, Dirk Hackenberger¹, Tara Sadoway¹, Ross MacPhee⁵, Duane Froese², Hendrik Poinar¹



University





1. Introduction

• Targeted enrichment is a powerful means of improving the fraction of informative DNA from high-throughput sequencing, and can be used to capture entire genomes from highly degraded ancient remains.



• PCR metabarcoding is still the predominant means of isolating environmental DNA for reconstructing modern and palaeo-ecosystems. But, this approach is limited by:

- compounding PCR biases,
- an inability to recover small fragments or authenticate ancient molecules,
- limited genetic target loci.

• Here, we demonstrate the effectiveness of **PalaeoChip** enrichments utilizing Yukon permafrost cores dated to the Pleistocene-to-Holocene transition¹⁻³.

• **PalaeoChip** baits are designed to simultaneously capture chloroplast barcoding loci from 2500+ northern plant species and whole mitochondrial genomes from hundreds of Holarctic animals.

• We observe a 14.6-fold increase in ecological informative plant and animal DNA with our cold spin extractions and PalaeoChip compared with a commercial soil extraction kit, and a 22.6-fold increase compared with PCR metabarcoding.



Yukon coring sites and landscape extent during Last Glacial Maximum (LGM, 26–19 kya).

2. Methods



- Exterior surface removed with sterilized chisel while frozen.
- **2.** Liquid nitrogen cooled hollow drill bit used to repeatedly subsample core section.

3. Subsamples homogenized by stirring.



Metagenomic comparison of Bear Creek permafrost core (30,000 cal-BP) (left). Libraries⁴⁻⁵ compared in MEGAN⁶⁻⁷ with absolute read counts using logarithmically scaled bubble charts indicating the number of unique reads assigned to each taxon node. MapDamage⁹ plots (*right*) show characteristic ancient DNA damage.

4. Discussion

σ

ad

- Our novel inhibitor removal procedure paired with Dabney⁸ purifications and the **PalaeoChip Arctic 1.0** targeted enrichment bait-set outperformed a sedaDNA commercial extraction kit across all extraction replicates, as well as outperforming a PCR metabarcoding approach and shotgun sequencing.
- With **PalaeoChip**, we can simultaneously capture a complex set of highly degraded environmental DNA from a range of informative loci (including across whole organelle genomes) while also evaluating ancient DNA authenticity⁹.
- There is still room for further optimization (*right*) of our environmental DNA inhibitor removal procedure to increase

Balancing inhibitor removal with DNA retention during sedaDNA extraction





PowerBead tubes with a custom digest solution used to lyse sedaDNA.

• Proteinase K added after vortexing, then incubated @ 35°C overnight while spinning.

• Lysis supernatant pipetted into high-volume binding buffer. • Centrifuged overnight @ 4°C. • Supernatant decanted, purified following Dabney et al.⁸





paired-end sequencing

its efficiency for other sedimentary contexts. We also intend on further optimizing the **PalaeoChip Arctic** baits with an

expanded set of ecologically informative target organisms,

and for other geographic regions.



Murchie, Tyler J., Melanie Kuch, Ana Duggan, Marissa L. Ledger, Kevin Roche, Jennifer Klunk, Emil Karpinski, Dirk Hackenberger, Tara Sadoway, Ross MacPhee, Duane Froese, Hendrik Poinar 2019. PalaeoChip Arctic1.0: An optimised eDNA targeted enrichment approach to reconstructing past environments. bioRxiv 730440.

to library adapt DNA inhibitors

> Inhibition (DNA dependent and independent)

5. Acknowledgements and References

Our thanks to **Alison Devault** at Arbor Biosciences for her help with designing the baits. Also thanks to everyone in Brian Golding's bioinformatics group, all members of the McMaster Ancient DNA Centre, and the McMaster Anthropology department for their continued support. Thanks to SSHRC, NSERC, the W. Garfield Weston Foundation, and ACUNS for their financial support.

1. T. R. Sadoway, A Metagenomic Analysis of Ancient Sedimentary DNA Across the Pleistocene-Holocene Transition, Master's thesis, McMaster University (2014).

2. M. E. Mahony, 50,000 years of paleoenvironmental change recorded in meteoric waters and coeval paleoecological and cryostratigraphic indicators from the Klondike goldfields, Yukon, Canada, Master's thesis, University of Alberta (2015).

3. V. M. D'Costa et al., Antibiotic resistance is ancient. Nature. 477, 457–461 (2011).

4. M. Meyer, M. Kircher, Illumina sequencing library preparation for highly multiplexed target capture and sequencing. *Cold Spring Harb. Protoc.* 5 (2010), doi:10.1101/pdb.prot5448.

5. M. Kircher, S. Sawyer, M. Meyer, Double indexing overcomes inaccuracies in multiplex sequencing on the Illumina platform. Nucleic Acids Res. 40, 1-8 (2012).

6. D. H. Huson et al., MEGAN Community Edition - Interactive Exploration and Analysis of Large-Scale Microbiome Sequencing Data. PLoS Comput. Biol. 12 (2016), doi:10.1371/journal.pcbi.1004957.

7. D. H. Huson, A. F. Auch, J. Qi, S. C. Schuster, MEGAN analysis of metagenomic data. Genome Res. 17, 377–386 (2007).

8. J. Dabney et al., Complete mitochondrial genome sequence of a Middle Pleistocene cave bear reconstructed from ultrashort DNA fragments. Proc. Natl. Acad. Sci. 110, 15758–15763 (2013).

9. H. Jónsson, A. Ginolhac, M. Schubert, P. L. F. Johnson, L. Orlando, MapDamage2.0: Fast approximate Bayesian estimates of ancient DNA damage parameters. Bioinformatics. 29, 1682-1684 (2013).